



# Brain Profiler

## BRAIN PROFILER USER GUIDE

<http://www.brainprofiler.com/>






# Brain Profiler

Brain Profiler is an app you share with your clinician. You download it after you agree with your clinician on a password that you will share, and allows you to use the application and keep in contact with your clinician

On your first visit your clinician will explain to you how the application works, thereafter consecutive follow-up on your treatment can be done with the help of the application making your follow-up and treatment more personalized and effective





1. Slide to rate

Touch and slide,  
positioning the  
cursor on the  
correct statement  
about your feeling



2. Marks pages  
progress 9 in all

After rating, the  
page turns  
automatically to the  
next page, you can  
also slide the screen  
to go back and forth  
between pages

Carrier 16:59

Brain Profiler

So how are you today?

Nothing interests,  
nothing is important,  
cannot do anything

Disinterested hard to  
act

Unmotivated

None



# Brain Profiler

Results of your report are available for your clinician to make personalized decisions about your treatment. He may reschedule you a visit, or contact you by visual (e.g., Skype) auditory (phone) or SMS media.



## IMPORTANT NOTICE

Brain Profiler outputs your clinical condition as diagnosed traditionally – however it also proposes a brain-related prediction of the disturbance. The brain-related proposal is literature-based but has not been validated yet and is in research phase.